

Spot test name	Cobalt	Marquis	Froehde's	Meckes	Dilli	Erlich's	Duquenois
Drug Name							
Acetaminophen	(-)		Light Blue	Faint yellow			
Amphetamine		Orange > brown	(-)	(-)			
Aspirin		raspberry after 10 mins	violet				
Barbituates					Purple		
Benactyzine		Orange > green > blue	Orange > brown				
Benzocaine	(-)	(-)	(-)	(-)	(-)		
Bufotenine			Purple				
Chlorcyclizine			faint yellow				
Cocaine	blue	wine red w/heat					
Codeine	blue	Red-violet > blue-violet	green > red-brown	green > blue-green			
Cyclizine			Yellow				
Demerol		Orange					
DET			(-)			Purple	
Dilaudid		Red-violet > blue-violet	Violet > blue				
Diphenhydramine		Yellow	Yellow	Yellow			
DMT			(-)			Purple	
DOE			Yellow > greenish yellow				
Dolophine (Methadone)	blue	(-)	Yellowish green > green	Pink			
DPT			faint violet				
fentanyl		faint orange					
Hashish							Red-violet
Heroin		Purple	Purple > green	blue-green			
Ibogaine		grey > orange	Brown > purple > blue				
Lidocaine	blue						
LSD		orange > brown > purple	Olive green > blue green			Purple	
Marijuana							red-violet
MDA		black	Brown > purple				
MDMA		purple > black	green > black	green > blue-green			
Mescaline		orange	Yellowish brown				
methamphetamine		bright orange	(-)				
Methylphenidate			(-)				
Methypylon			light blue				
Morphine		Purple	Purple > slate				
MPP		faint orange					
MPPP		red					
MPTP		red					
Numorphan		Purple > light green	blue > green				

Opium		Purple	(-)				
Pantopon		Purple					
Papaverine		pink > light blue					
Phencyclidine			(-)				
Phenmetrazine			(-)				
procaine	blue						
Propoxyphene		black	brown > purple				
Psilocin			green > blue-green				
Psilocybin		orange	Olive green > yellow				
STP		Yellow > brown	Yellow > greenish yellow				
TMA			Blue & brown streaks > brown				